PARI PEAK FLOW AIR

ERWACHSENE & KINDER

CE
0086

Peak Flow Meter

Ref. No. 43904

EN ISO 23747

Model 4300

5 099169 439049

PZN-11024966

Manufacturer: Vitalograph, Ennis, Ireland

ENG Use, Cleaning & Regulatory Notices

Setting The Action Zones

The zone markers may be set with a stylus by the healthcare professional to give the personal boundaries for implementing the asthma action plan.

How To Use The Peak Flow Meter

- Slide the pointer to the bottom of the scale
- Hold the device horizontally in front of your mouth
- Inhale as deeply as possible
- Seal your lips firmly around the mouthpiece
- Blow out for a second or longer as HARD and as FAST as possible
- Take care that the tongue or teeth do not obstruct the airflow. A 'spitting' action or coughing will give false readings
- Take note of the reading and slide the pointer back to the bottom of the scale
- Repeat this three times and record the highest reading
- Take action according to your personal asthma management plan

Cleaning

The outer surfaces should be cleaned every week, more often if necessary. The use of an ordinary alcohol wipe (IPA 70%) is recommended, with special attention to the mouthpiece area.

Indications For Use/Intended Use

Measure peak expiratory flow (PEF), monitor progress in treatment of asthma, determine the severity of the asthma, check response to treatment during an acute asthma episode and detect worsening lung function and avoid a possible severe asthma flare-up.

Use In Clinic

When used for multiple subjects fit the mouthpiece adaptor (available as an accessory) and always use a disposable SafeTway® mouthpiece. The Peak Flow Meter should be externally cleaned between each patient. Medical facilities may use a spirometer in series to regularly check the accuracy of this device.

Warnings, Advisory Notices/Contraindications

Note: Please read all this information before using this Peak Flow Meter.

- This device is intended to help you assess how well your asthma is under control, by measuring your peak flow which can then be monitored over time
- Never attempt to dismantle the unit. This can cause faulty Peak Expiratory Flow scores
- If used at home symptoms must take precedence over device measurements
- If you think that the device is not reading correctly, advise your healthcare professional immediately
- If the device is used for longer than its specified life, the accuracy of the device may deteriorate
- Store in a clean dry place

Warranty

Your device is guaranteed for one year*. Replace if it is faulty, otherwise replace the unit every three years.

* Excepting accidental/transit damage or inappropriate use of the device.

CE Notice

Marking by the  symbol indicates compliance of this device to the Medical Devices Directive of the European Community. Such marking is indicative that the device meets or exceeds the referenced technical standards.

FDA Notice

Caution: Federal Law restricts this device to sale by, or on the order of a physician.

Declaration Of Conformity

Product: Vitalograph Model 4300
Vitalograph hereby ensures and declares that the above product associated with this user manual, is designed and manufactured in accordance with the following QMS regulations and standards:

- European Medical Devices Directive 93/42/EEC, as amended. This device is classified as 1m per Annex V per Annex VII section 5 of the MDD
- Canadian Medical Device Regulation {CMDR SOR/98-282}
- FDA Quality System Regulation (QSR) 21 CFR 820;
- EN ISO 13485. Medical devices. Quality management systems. Requirements for regulatory purposes

0086

Certifying Body (For 93/42/EEC and CMDR):
British Standards Institute, Certificate Nos.
CE 00772, CE 85553, MD 82182.

Technical Specifications

Size: 135.0mm x 54.2mm

Weight: 78.4g (Net)

Package Contents: Peak Flow Meter x 1, User Instruction Leaflet x 1

Product: Peak Flow Meter

Model Number: 43904

Material: ABS plastic

Accuracy:

± 10 L/min or ± 20% of the reading

Repeatability:

± 10 L/min or ± 5% of the reading

Highest resistance to flow:

0.00377 kPa/L/min @ 720 L/min

Measurement Range: 50-800 L/min BTPS

Operating Temperature Range:

ATS/ERS: Limits: 17-37°C.

Design Limits: 10-40°C

Storage Conditions:

Temperature: 0-50°C

Relative Humidity: 10%-95%

Peak Flow Meter performance standards:

EN ISO 23747:2015; ATS/ERS:2005;

AS/NZS4237


Frequency: Profile A/B difference less than

15 L/min/12% (Annex B, EN ISO 23747:2015)

Designed & Manufactured Under:

EN ISO 13485, FDA 21CFR820, CMDR & JPAL

Signed on behalf of Vitalograph (Ireland) Ltd.

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Gezeichnet im Namen von Vitalograph (Ireland) Ltd.