

## TREATMENT - SQ SQUARE TROLLEY Q638 B

**DIMENSIONS** 480 x 480 mm

**ADEQUATE HEIGHT, ERGONOMICS** 980 mm (work surface – floor).

### WORK SURFACE

- Polyethylene, PE. Asepsis, hygiene, easy-to-clean.
- Withstands alcohol, detergents, disinfectants, chlorhexidine.
- Protection on all 4 sides. Pleasant to the touch, rounded shapes.

### MONOBLOC STRUCTURE

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

### SQ SQUARE DRAWERS

- Polypropylene, PP. 2 colours - B Blue, W White.
- Total opening.
- Steel ball-bearing slides with self-closing system.
- Front stop, it protects the slide against strong impacts when closing.
- Load capacity per drawer: 10 kg.
- Ergonomic handle.
- Drawers with white cards, removable.

### EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION

- Elevated push handle. Ergonomics, height 1020 mm.
- Double wheels Ø 100 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at push handle side.
- Wheels Ø 125 mm ref. AD12 (optional).



## Q638 B

SQ065 B.AST.32R.46.50

### ACCESSORIES

AST ELEVATED PUSH HANDLE

32R SIDE SUPPORT

46 3 GLOVES BOX HOLDER

50 13 LITRES BIN