

STORAGE FURNITURE - CATHETER STORAGE TROLLEY 706 DN.60CAT

DIMENSIONS 640 x 480 mm

HEIGHT 1610 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

5 EXTRACTABLE SUPPORTS

- Adjustable ISO shelf with 5 extractable supports for catheters (25 hooks).
- 5 identification cards ISO norm 7810, international standard.

ISO PANELS

- Flexible organisation inside the trolley.

2 DOORS

- Total opening 270°.
- Phenolic, 10 mm thickness.
- Lock with 4 keys, multiple (2 trolleys = 2 equal locks).
- Card holder with white cards, removable, left door.
- DN Grey.
- DB Blue.

STATIC TROLLEYS

- Double wheels Ø 100 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at the front part.



706 DN.60CAT

ACCESSORIES

60CAT ADJUSTABLE ISO SHELF FOR CATHETERS