



VELA Medical

SKAI APOLLO



GREY - 057



BLACK - 052

Top coat: 88,7% PVC :: Backing 11,3% polyester ::
50.000 cycles DIN53359 :: EN 1021-1 :: EN 1021-2

SKAI AMBLA



TINTO BLUE - 023

Top coat: 100% PVC :: Backing 100% cotton ::
300.000 cycles Martindale :: EN 1021-1 :: EN 1021-2

strack Patienten-
und Pflegebedarf

Strack AG
Ebnatstrasse 125
8200 Schaffhausen

052 630 31 60
www.strack.ch
info@strack.ch