

360-0300
SQM EMERGENCY TROLLEY FUNCTIONALITY
FUNCIONALIDAD CARRO DE EMERGENCIA SQM

Square extensive mobility trolleys for use in hospitalization units.
Dimensions 480 x 480 mm.
Adequate height, ergonomics: 1005 mm (work surface - floor).
Square drawers. Heights 57, 130 mm.
Manual locking for drawers with disposable seals.
This safety system guarantees that the content of the trolley is intact at the critical time of use.

Carros de alta movilidad cuadrados para uso en unidades de hospitalización.
Medidas 480 x 480 mm.
Altura adecuada, ergonomía: 1005 mm (superficie de trabajo - suelo).
Cajones cuadrados. Alturas 57, 130 mm.
Bloqueo manual de cajones con precintos desechables. Este sistema de seguridad garantiza que el contenido del carro se encuentra intacto en el momento crítico de utilización

SERUM
SUEROS

CARDIAC
MESSAGE BOARD
TABLA DE MASAJE
CARDÍACO

DEFIBRILLATOR
2 ROTATION AXLES
DESFIBRILADOR
2 EJES DE GIRO

ELEVATED PUSH
HANDLE
TIRADOR ELEVADO

PROTECTION
ON ALL 4 SIDES
PROTECCIÓN 4 LADOS

IDENTIFICATION
IDENTIFICACIÓN

DISPOSABLE
CONTAINER

GAS BOTTLE
Ø 120 / 140
BOTELLA DE GAS
Ø 120 / 140


CONTENEDOR
DESECHABLE

ERGONOMIC
HANDLE
ASA
ERGONÓMICA

SEAL LOCKING
BLOQUEO PARA PRECINTOS

ERGONOMIC PEDAL
PEDAL ERGONÓMICO

2 SQUARE DRAWERS 57 mm
(Optional)
2 CAJONES CUADRADOS 57 mm
(Opcional)



Q260 R

6 DMQ Per drawer Por cajón
16 compartments compartimentos