

**TRANSPORT TROLLEYS - 500 SERIES
566**

DIMENSIONS 900 x 630 mm
HEIGHT 1310 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.
- Interior shelves for instrument boxes of 600 x 300 mm.

2 DOORS

- Stainless steel AISI-304.
- Total opening 270°.
- Lock with 4 keys.

LOAD CAPACITY

- 150 kg.

STERILISATION

- Non-airtight doors, the interior of the trolley is not isolated from contact with the outside air.

EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION

- Elevated push handle.
- Dual lock wheels Ø 125 mm. Soft tread, silent.
- 2 antistatic (electrical conductive). 2 with brake at the right side of the trolley, opposite the handles.
Pedal position: Superior - fixed wheel, easy driving.
Central - swivel wheel.
Inferior - locked wheel.

