

**TRANSPORT TROLLEYS - 500 SERIES**  
**536**

**DIMENSIONS** 900 x 630 mm

**HEIGHT** 930 mm

**MONOBLOC STRUCTURE**

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

**PROTECTION**

- 3 protection rods for each self.

**LOAD CAPACITY**

- 3 shelves, 150 kg.

**EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION**

- Elevated push handle.
- Double wheels Ø 125 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at push handle side.

