

TRANSPORT TROLLEYS - 500 SERIES
543

DIMENSIONS 900 x 500 mm

HEIGHT 1280 mm

MONOBLOC STRUCTURE

- Stainless steel AISI-304. Fine satin finish, asepsis, hygiene.

PROTECTION

- 3 protection rods for each self.

LOAD CAPACITY

- 4 shelves, 150 kg.

EXTENSIVE MOBILITY TROLLEYS, DAILY CONDUCTION

- Elevated push handle.
- Double wheels Ø 125 mm. Soft tread, silent.
- 2 antistatic (electrical conductive), 2 with brake at push handle side.

