


ADMIRAL GRAPH CHART

MINIMUM ANCHOR SETBACK = 14.5"

MAXIMUM ANCHOR SETBACK = 19"

DRAW YOUR PROFILE ON THE GRAPH USING THE FRONT ANCHOR LOCATIONS AS YOUR START POINT. ALL GRAPH LINES ARE IN 1" INCREMENTS.



strack Liftsysteme

Strack AG
Ebnatstrasse 125
8200 Schaffhausen

Tel. 0800 600 500
Fax 0800 600 900

www.strack.ch
info@strack.ch

ZSR-Nr. K 7514.14
NIF-Nr. 366 824
MWST-Nr. 567 970

PRODUCT:		ADMIRAL	
PART #:		F-ADMRL	
SCALE:	1:9	DATE:	7/15/2019
DRAWN BY:	KR	REVISION:	03