



- FEATURES**
- ◆ UL CERTIFIED ◆ 450 LB. (204 KG.) WEIGHT CAPACITY ◆ SETBACK RANGE 14 1/2" TO 19" (37 - 48 cm) ◆ FIELD REVERSIBLE ◆ ADJUSTABLE BASE PLATE
  - ◆ FLIP UP ARMREST ◆ ALL STAINLESS STEEL CONSTRUCTION ◆ DURABLE POWDER COAT FINISH ◆ RECHARGEABLE BATTERY & WALL MOUNT CHARGER
  - ◆ FIVE YEAR STRUCTURAL WARRANTY ◆ FIVE YEAR PRO-RATED ELECTRONICS WARRANTY

CLICK HERE TO VIEW THIS PRODUCT:  
[HTTPS://AQUACREEKPRODUCTS.COM/ADMIRAL-POOL-LIFT/](https://aquacreekproducts.com/ADMIRAL-POOL-LIFT/)

**ADMIRAL**  
U.S. PATENT NO. 7,249,386 B2 & D507,769 S  
PART NO. F-ADMRL

MATERIAL: 304 STAINLESS STEEL  
FINISH: POWDER COAT - WHITE  
APPROX. WEIGHT: 119.86 LBS.

**TOLERANCES**  
DIMENSIONS ARE IN INCHES  
UNLESS OTHERWISE SHOWN  
FRACTIONAL: ..... ±1/16  
ANGULAR: ..... ±1°  
TWO PLACE DECIMAL: .... ±.01  
THREE PLACE DECIMAL: ... ±.005



9889 GARRYMORE LANE  
MISSOULA, MT 59808  
(406) 549-0769  
FAX (406) 549-2602  
www.AquaCreek.com

**PROPRIETY AND CONFIDENTIAL INFORMATION**  
THIS DRAWING AND ALL INFORMATION CONTAINED HEREIN IS THE SOLE PROPERTY OF AQUA CREEK PRODUCTS LLC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN PERMISSION FROM AQUA CREEK IS PROHIBITED.  
\*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
P1	INITIAL RELEASE	4/5/2021	J.MAXWELL

**PRODUCT NAME:** Admiral™ Pool/Spa Access Lift

**CSI SECTION:** 13 11 46 Swimming Pool Accessories

**APPLICATIONS:**

- In Ground Pools & Spas with flush pool decks having up to **2-inch (5.08cm)** high curbs.

**SUBMITTALS:**

- Alternate make or model lifts will not be considered unless equal to the specified product in every respect. Alternative products must be submitted for approval no less than ten (10) calendar days prior to the bid due date. Submittal data must include complete documentation relating to all the specified features and include all manufacturer's sales literature, specification sheets, installation/operation manuals, engineered drawings and documents that support the strength, service capacity, and the ability to meet the appropriate Federal Regulatory Standards associated with the product.

**QUALITY ASSURANCE:**

- The product must be a standard product of a company registered with the United States Food and Drug Administration (FDA) as a class 1 medical device manufacturer and that is regularly engaged in the manufacture of commercial grade swimming pool and spa accessibility equipment.
- The product shall be structurally capable of supporting a **450 pound (204.12kg)** load capacity and static load tested to 1.5 times the rated load.
- The product shall be certified to be compliant with UL test standard 60335-2-1000 Standard of Household and Similar Appliances: Particularly Requirements for Electrically Powered Pool Lifts.
- The product shall be tested and confirmed to comply with design standards outlined in the United States ADA Accessibility Standards for Accessible Designs.

**SPECIFICATIONS:**

- The products main structure shall be constructed of type 304 stainless steel; no dissimilar metallic components shall be allowed. The product shall be finished with a minimum 4-mil thick, UV resistant powder coat finish. All fastening hardware shall be of type 316 stainless steel. Product shall be supplied with a Linak Actuator suitable for medical use, a separate wall mount battery charger, a hand-held controller, and a 24-volt rechargeable battery.
- The chair of the product shall be of a rigid one-piece, polyethylene, molded design. The minimum seat width shall be **18-inches (45.72cm)** and the seat back shall be a minimum height of **17.5-inches (44.45cm)**. The seat shall feature dual flip-up arms and supplied with a 2-inch (5.08cm) wide adjustable safety belt with hook and loop fasteners.
- The seat height of the product will be adjustable from a load height of **18-1/4-inches (46.36cm)** to **24-1/4-inches (61.60cm)** in **2-inch (5.08cm)** increments to accommodate different height wheelchairs.
- The product shall be capable of accommodating a deck to water distance (water draft) of up to **8-inches (20.32cm)** and capable of an anchor set back distance of up to **19-inches (48.26cm)**.
- Product can be field reversible and includes an adjustable and removable footrest that moves with the seat.

**WARRANTY:**

- Product shall have a minimum 5 year structural and 5 year pro-rated electronic warranty.

**LIFT PART NUMBERS:**

Item No	Order No	Description
F-ADMRL	INV01513	Admiral Lift White Powder Coat w/Blue Seat
F-ADMRL-T	INV01757	Admiral Lift White Powder Coat w/Tan Seat
F-ADMRL-G	INV01756	Admiral Lift White Powder Coat w/Gray Seat
F-ADMRL-W	INV01755	Admiral Lift White Powder Coat w/White Seat
F-ADMRL-C	INV01907	Admiral Lift (Specify Seat and Powder Coat Colors)

**ANCHORING OPTIONS:**

Item No	Order No	Description
F-04CAJP	INV01411	Admiral Lift Anchor for Standard Concrete Decks
F-08CAJP	INV01414	Admiral Lift Anchor for Paver Applications
F-046QAB-AT1	INV01406	Admiral Lift Anchor for Standard Concrete Decks with Quick Attach Knobs
F-40220-AT1	INV01444	Admiral Lift Anchor for Wood Deck Applications

**PRODUCT OPTIONS:**

Item No	Order No	Description
F-0904CVR	INV01990	Admiral Lift Protective Seat Cover in Blue
F-120PPC	INV00525	Admiral Lift Protective Lift Cover in Blue
F-120PPC-T	INV00528	Admiral Lift Protective Lift Cover in Tan
F-120PPC-G	INV00527	Admiral Lift Protective Lift Cover in Gray
F-044SCH	INV01693	Admiral Lift Solar Charger
F-422PLH-2	INV01455	Admiral Lift Adjustable Head Rest
F-423CS	INV01457	Admiral Lift Chest Strap
F-428HT	INV01458	Admiral Lift Transport Cart
F-LRB	INV01422	Admiral Lift Pull Out Leg Rest Blue
F-LRT	INV01952	Admiral Lift Pull Out Leg Rest Tan
F-LRG	INV01954	Admiral Lift Pull Out Leg Rest Gray
F-LRW	INV01953	Admiral Lift Pull Out Leg Rest White
F-019CA	INV01379	Admiral Lift Cycle Attachment

Aqua Creek Products LLC  
9889 Garrymore Lane  
Missoula, MT 59808

[www.aquacreekproducts.com](http://www.aquacreekproducts.com)

Tel: 888-687-3552 / 406-549-0769

Fax: 406-541-8958



9889 GARRYMORE LANE  
MISSOULA, MT 59808  
(406) 549-0769  
FAX (406) 549-2602  
www.AquaCreek.com

**PROPRIETY AND CONFIDENTIAL INFORMATION**

THIS DRAWING AND ALL INFORMATION CONTAINED HEREIN IS THE SOLE PROPERTY OF AQUA CREEK PRODUCTS LLC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT WRITTEN PERMISSION FROM AQUA CREEK IS PROHIBITED.

\*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

© 2020 AQUA CREEK PRODUCTS. ALL RIGHTS RESERVED.

MODEL  
**P1**  
REVISION

**ADMIRAL - F-ADMRL**

51-163-220

SCALE  
1:14

SHEET  
**2**