


EZ & POWER EZ 10" ANCHOR




strack Liftsysteme

Strack AG
Ebnatstrasse 125
8200 Schaffhausen


Tel. 0800 600 500
Fax 0800 600 900

www.strack.ch
info@strack.ch

ZSR-Nr. K 7514.14
NIF-Nr. 366 824
MWST-Nr. 567 970



SECTION A-A



GROUNDING LUG

SPECIFICATIONS:

RATED CAPACITY: 400LBS [181kg]

MAXIMUM SETBACK: 20" [508mm]

PRODUCT:		EZ 10" ANCHOR	
PART #:		F-30220	
SCALE:	1:3	DATE:	6/12/2017
DRAWN BY:	KR	REVISION:	A

