

POWER EZ 2

EZ SERIES

ITEM # F-PEZPL2



ABOUT

With all the ease and simplicity of the EZ 2 Lift, the Power EZ 2 is an even better solution for affordable and compact lifts! With a manual 360 degree rotation, this lift utilizes a battery powered actuator to raise and lower users with confidence. In addition to a setback range of 12"-20", the adjustable sling seat allows for 20" of wall clearance.

Best of all, this lift is easily portable for storing and relocation! All this, plus the best warranty in the industry, makes our EZ line an easy choice for your accessibility needs.

FEATURES:

- 400 lb (181 kg) weight capacity
- Durable mesh sling
- Adjustable stainless steel chains
- Stainless steel construction with powder coat finish
- Limited five (5) year warranty
- Ships small parcel, some assembly required
- Made in USA

ACCENT COLORS



EZ Series does not come standard with color options.

strack Liftsysteme

Strack AG
Ebnatstrasse 125
8200 Schaffhausen

Tel. 0800 600 500
Fax 0800 600 900

www.strack.ch
info@strack.ch

ZSR-Nr. K 7514.14
NIF-Nr. 366 824
MWST-Nr. 567 970

Affordable recreation equipment for people of all abilities.

1.888.687.3552 | aquacreek.com | sales@aquacreek.com

POWER EZ 2

APPLICATIONS

This lift works on in ground and partially raised pools and spas. To see if the EZ 2 will work for you, visit our website and download the application guide, or contact our friendly sales team for assistance. This lift can also be anchored to pool decks made of concrete or wood. Choose the appropriate anchor system for your needs. Anchors sold separately.


Standard Anchor



Surface Mount Anchor



Wood Deck Anchor



For access to CADs, Revit drawings and more visit our Partner Login page on our website and become an Aqua Creek partner.

OPTIONS / ACCESSORIES




- Anchors for concrete, or wood deck applications
- Hard seat option
- Powder coat and accent color options to match almost any decor
- Hard seat cover
- Cover to protect your investment from the elements*

* Color and material options available

PREMIUM COVER



SOLAR CHARGER



Accessorize your lift with a cover to keep it looking like new and a solar charging station for an environmentally friendly power option.



AVAILABLE SEAT OPTIONS

Seat options in both hard and soft varieties to suit any needs.

