

MIGHTY 400



Item# F-MTY400
400 lb weight capacity
Up to 12" deck to water

The new Mighty Lifts are UL Certified, ADA compliant lifts designed to provide versatility and accessibility solutions with greater water drafts and setbacks up to 44 inches.

The Mighty series features an adjustable height seat pole allowing 4" of additional adjustment in 2" increments.

In addition the Mightys' can rotate 360 degrees in either direction. Designed for in ground anchoring, the Mighty series will accommodate existing "standard" Scout 2 and Revolution concrete anchors. Created for its adaptability, the Mightys' can be purchased with a pedestal for additional wall clearance or unusual gutter profiles. Pedestals come in 1', 2', 3', and 4' options.

The Mighty series is ideal for easy and affordable ADA compliance. US Patent Pending.

All lifts come standard in white with blue seat and footrest. Additional color options are available for both plastics and powder coat.

FEATURES:

- UL certified
- 400 LB weight capacity
- Up to 12" deck to water
- 360 degree rotation
- Setback up to 44"
- Adjustable seat height
- Direct drive rotation
- Flip-up armrests
- Stainless steel construction
- Durable powder coat finish
- Five year pro-rated electronics warranty
- Five year structural warranty

pedestal accessory



U.S. Patent No.
D507,769S

ACCENT COLORS

Available in a variety of accent and powder coat colors



*Blue seat standard (not a powder coat option)

OPTIONS/ACCESSORIES:

F-808SA	Anchor kit, standard 1.9" ID x 6" deep
F-808SA-10	Anchor, 10" for pavers
F-77SSA	Anchor, Semi-Recessed, Shallow Concrete
F-MTLC	Cover, standard, blue
F-MTLC-T	Cover, standard, tan
F-MTLC-G	Cover, standard, gray
F-422PLH-2	Headrest
F-LRB	Pull-out leg rest assembly (blue, tan, gray, or white)
F-048SCH	Solar charging station
F-MTTC	Mighty series transport cart
F-MTY1PD	Pedestal, 1 foot tall
F-MTY2PD	Pedestal, 2 foot tall
F-MTY3PD	Pedestal, 3 foot tall
F-MTY4PD	Pedestal, 4 foot tall

PROUDLY SOLD BY:

strack Liftsysteme

Strack AG
Ebnatstrasse 125
8200 Schaffhausen


Tel. 0800 600 500
Fax 0800 600 900

www.strack.ch
info@strack.ch

ZSR-Nr. K 7514.14
NIF-Nr. 366 824
MWST-Nr. 567 970

THE MIGHTY 400 POOL LIFT


PRODUCT SPECIFICATIONS




APPLICABLE POOL PROFILES

WITHOUT PEDESTAL


ROLL OUT




PARTIALLY RECESSED



FULLY RECESSED




BULL NOSE




WITH PEDESTAL


RAISED POOL / SPA



FLORIDA ROLL OUT



PARTIALLY RECESSED
W/ CURB



To confirm specific application go to www.aquacreek.com for our Pool Lift Selector
Become a member of our new Dealer and Architect page. Access to CAD blocks,
CSI specs and sales materials.