


# PORTABLE PRO POOL 2 GRAPH CHART

FRONT EDGE OF THE CART CAN BE PLACED A MAXIMUM OF 6" FROM THE POOL'S EDGE.

DRAW YOUR PROFILE ON THE GRAPH USING THE FRONT ANCHOR LOCATIONS AS YOUR START POINT. ALL GRAPH LINES ARE IN 1" INCREMENTS.



**strack** Liftsysteme

Strack AG  
Ebnatstrasse 125  
8200 Schaffhausen

Tel. 0800 600 500  
Fax 0800 600 900

www.strack.ch  
info@strack.ch

ZSR-Nr. K 7514.14  
NIF-Nr. 366 824  
MWST-Nr. 567 970

  
**Aqua Creek Products**

PRODUCT: <b>PORTABLE PRO POOL 2</b>	
PART #: <b>F-PPP2</b>	
SCALE: <b>1:12</b>	DATE: <b>7/24/2020</b>
DRAWN BY: <b>KR</b>	REVISION: <b>01</b>