

REVOLUTION

REVOLUTION SERIES

ITEM # F-702RLNA



U.S. Patent No. D507,769S



ABOUT

Thanks to its incredible diversity, The Revolution Pool Lift is making waves in the world of aquatics accessibility. From its 500 pound weight capacity to its adjustable seat pole and submersible 4-button handset, it's no wonder how this user-friendly lift got its name. While the whisper quiet 360 degree rotation guides you safely to and from the water, you can worry more about enjoying your day at the pool, and less about how you're going to get there.

Found in hundreds of pools across the globe, The Revolution Lift is designed with the same stainless steel quality that you can trust from all Aqua Creek products. Ask our friendly sales team how The Revolution can elevate your experience today!

FEATURES: ⚙️

- ADA Compliant & UL Safety Certified
- 500* lb (226 kg) weight capacity
- 360 degree rotation
- Adjustable leg rest
- Wider, more comfortable 22" seat
- Stainless steel construction with powder coat finish
- Limited five (5) year warranty
- Ships common carrier, some assembly required
- Made in USA

*400 lb capacity for deep draft version

*350 lb (158 kg) capacity with gurney, spine board, and wheelchair attachments

ACCENT COLORS 🎨

Available in a variety of accent and powder coat colors



*Blue seat standard (not a powder coat option)

strack Liftsysteme

Strack AG

Ebnatstrasse 125

8200 Schaffhausen

Tel. 0800 600 500

Fax 0800 600 900

www.strack.ch

info@strack.ch

ZSR-Nr. K 7514.14

NIF-Nr. 366 824

MWST-Nr. 567 970

Affordable recreation equipment for people of all abilities.

1.888.687.3552 | aquacreek.com | sales@aquacreek.com

REVOLUTION


APPLICATIONS

This lift works on in ground and some partially raised pools and spas. To see if this lift works with your pool, visit our website and download the application guide for this product, or contact our friendly sales team. This lift can be anchored to pool decks made of concrete. Choose the appropriate anchor system for your needs. Anchors sold separately.

Standard Anchor



Semi-Recessed Anchor



For access to CADs, Revit drawings and more visit our Partner Login page on our website and become an Aqua Creek partner.

OPTIONS / ACCESSORIES



- Anchoring options for concrete and paver applications
- Chest strap for additional upper body support
- Cover to protect your investment from the elements*
- Cycle attachment to convert your lift into an aquatic bike
- Headrest for additional support
- Pull out leg rest for additional wall clearances
- Rescue spine board option
- Sling seat option
- Solar charging station to conserve energy
- Transport cart for easy removal from pool side for storage**

*Color and material options available


**Transport cart not usable with semi recessed anchor

PREMIUM COVER



Accessorize your lift with a cover to keep it looking like new and a solar charging station for an environmentally friendly power option.

SOLAR CHARGER



CUSTOM COLOR OPTIONS

Match your decor with a variety of custom color options.

