



SCOUT 2 POOL LIFT

U.S. PATENT NO. D507,769 S

PRODUCT WEIGHT: 210lbs [95kg]
 MAX WEIGHT CAPACITY: 375lbs [170kg]
 MIN ANCHOR SETBACK: 19"
 MAX ANCHOR SETBACK: 35"
 MAX WATER DRAFT: *10"
 SEAT HEIGHT ADJUSTMENT: 16" IN 4" INCREMENTS
 VERTICAL TRAVEL DISTANCE: 46"



PRODUCT:	SCOUT 2 LIFT		
PART #:	F-802SC2		
SCALE:	1:32	DATE:	8/15/2019
DRAWN BY:	KR	REVISION:	03



strack Liftsysteme