

SPA LIFT ELITE 68

SPA SERIES

ITEM # F-006SLE



U.S. Patent No. D507,769S



ABOUT

Crushing all limitations is the Spa Lift Elite, bringing you a lift specially designed for providing access to raised pools and spas, clearing walls from 21 to 68 inches high! With a 360 degree continuous rotation, adjustable seat height, and easy-to-use push button controls, this lift will completely change the way you look at pool and spa accessibility.

Ideal for swim spas and raised pools, this durable stainless steel design is backed by the industry's best warranty. The Spa Lift Elite comes in a number of attractive accent colors that will match almost any decor, and more customized colors are also available upon request!

FEATURES:

- ADA compliant & UL Safety Certified
- 300 lb (136 kg) weight capacity
- 360 degree rotation
- Adjustable leg rest & seat height
- Includes 24v rechargeable battery and charger
- Stainless steel construction with powder coat finish
- Limited five (5) year warranty
- Ships common carrier, some assembly required
- Made in USA

ACCENT COLORS

Available in a variety of accent and powder coat colors



*Blue seat standard (not a powder coat option)

strack Liftsysteme

Strack AG
Ebnatstrasse 125
8200 Schaffhausen

Tel. 0800 600 500
Fax 0800 600 900

www.strack.ch
info@strack.ch

ZSR-Nr. K 7514.14
NIF-Nr. 366 824
MWST-Nr. 567 970

Affordable recreation equipment for people of all abilities.


1.888.687.3552 | aquacreek.com | sales@aquacreek.com

SPA LIFT ELITE 68


APPLICATIONS

This lift works on raised pools and spas. To see if the Spa Lift Elite works on your pool, visit our website and download the application guide, or contact our friendly sales team. This lift can be anchored to pool decks made of concrete and wood. Choose the appropriate anchor for your needs. Anchors sold separately.

Concrete Anchor



Wood Deck Anchor



For access to CADs, Revit drawings and more visit our Partner Login page on our website and become an Aqua Creek partner.

OPTIONS/ACCESSORIES


- Anchors for concrete or wood deck applications
- Cycle attachment to convert your lift into an aquatic bike
- Cover to protect your investment from the elements*
- Headrest for comfortable head support
- Solar charging station to conserve energy
- Upper chest strap for added torso stability

* Color and material options available

PREMIUM COVER



SOLAR CHARGER



Accessorize your lift with a cover to keep it looking like new and a solar charging station for an environmentally friendly power option.



CUSTOM COLOR OPTIONS

Match your decor with a variety of custom color options.