




# SPA LIFT ULTRA


U.S. PATENT NO. D507,769 S



**PRODUCT WEIGHT: 180lbs [81kg]**  
**MAX WEIGHT CAPACITY: 400lbs [181kg]**  
**MIN ANCHOR SETBACK: 10"**  
**MIN WALL HEIGHT: 21"**  
**WALL THICKNESS @ MIN HEIGHT: 12"**  
**MAX WALL HEIGHT: 51"**  
**WALL THICKNESS @ MAX HEIGHT: 14"**  
**VERTICAL TRAVEL DISTANCE: 41 1/8"**



**LOADING POSITION**



**ABOVE WALL**



**strack** Liftsysteme

Strack AG      Tel. 0800 600 500      www.strack.ch      ZSR-Nr. K 7514.14  
 Ebnatstrasse 125      Fax 0800 600 900      info@strack.ch      NIF-Nr. 366 824  
 8200 Schaffhausen                     MWST-Nr. 567 970

PRODUCT: SPA LIFT ULTRA	
PART #: F-005SLU	
SCALE: 1:24	DATE: 8/22/2019
DRAWN BY: KR	REVISION: 04