SUPER POWER EZ

F - 1 9 4 P H B L - N A





ABOUT

Accessibility shouldn't be a burden - it should be easy! The Super Power EZ is here to provide you with a clean, compact design that is user friendly without compromising strength and durability. With all the simplicity of the EZ line, this lift has the ability to conquer greater obstacles with up to 24" setbacks and an impressive 58 inch wall clearance.

With a durable stainless steel design, rechargeable battery, and a number of accessory options, the Super Power EZ is an excellent option for affordable accessibility at home.

FEATURES: 🏠

- 400 lb (181 kg) weight capacity
- · Adjustable stainless steel chains
- Durable mesh sling
- · Rechargeable battery operation with charger
- Stainless steel construction with powder coat finish
- Limited five (5) year warranty
- Ships common carrier, some assembly required
- Made in USA

ACCENT COLORS





EZ Series does not come standard with color options.

strack Liftsysteme

Ebnatstrasse 125 8200 Schaffhausen

Fax 0800 600 900

info@strack.ch

www.strack.ch ZSR-Nr. K 7514.14 NIF-Nr. 366 824 MWST-Nr. 567 970



SUPER POWER EZ

APPLICATIONS

This lift works on raised pools and spas. To see if the Super Power EZ 2 will work for you, visit our website and download the application guide, or contact our friendly sales team for assistance.

This lift can be anchored to pool decks made of concrete only. Anchors sold separately.

Standard Concrete Anchor



For access to CADs, Revit drawings and more visit our Partner Login page on our website and become an Aqua Creek partner.

OPTIONS/ACCESSORIES



- Anchoring options for concrete deck applications only
- Cover to protect your investment from the elements*
- Four-point hanger bar option for added safety
- Hard seat option
- Various sling styles and sizes to choose from
- * Color and material options available

PREMIUM COVER







Accessorize your lift with a cover to keep it looking like new and a solar charging station for an environmentally friendly power option.



AVAILABLE SEAT OPTIONS

Seat options in both hard and soft varieties to suit any needs.

















